




| | |
|---|--|
| Name of Item | Advanced 7530 Wheelchair |
| Main target users' type of disabilities | Physical Impairment Chronic Medical Illness |
| Specification | Weight: 11 kg Load: 360 lb Seat width: 18" |
| Description | Foldable transit wheelchair |
| Product Information Page | N/A |



| | |
|---|--|
| Name of Item | 3M™ Adjustable Foot Rest FR430CB, Charcoal Black |
| Main target users' type of disabilities | Physical Impairment Chronic Medical Illness |
| Specification | Weight: 11 kg Load: 360 lb Seat width: 18" |
| Description | Footrest with adjustable height and inclination for more comfort while sitting |
| Product Information Page | Link |




| | |
|---|--|
| Name of Item | Standing Desks (1.0m / 1.2m) |
| Main target users' type of disabilities | Physical Impairment |
| Specification | <p>1.0m Desktop Size: 1000 x 600mm Height Range: 720mm – 1200mm</p> <p>1.2m Desktop Size: 1200 x 700mm Height Range: 645mm – 1295mm</p> |
| Description | Electric-powered, height adjustable tables to allow working while standing or to accommodate a wheelchair |
| Product Information Page | <p>1.0m: A6-N2-XS standing desk Free Delivery / 2 years Onsite Warranty (freemax.com.hk)</p> <p>1.2m: A6-N3 Standing Desk Free Delivery / 2 years Onsite Warranty (freemax.com.hk)</p> |



| | |
|---|--|
| Name of Item | MicroMLxS Receiver for BTE Hearing Aids AS10 Audioshoe for Naida UP BTE Hearing Aids |
| Main target users' type of disabilities | Hearing Impairment |
| Description | A Receiver device compatible with most behind-the-ear hearing instruments and Cochlear implants speech processors. To be used together with various FM Transmitters. |
| Product Information Page | Link |



| | |
|---|--|
| Name of Item | Roger MyLink |
| Main target users' type of disabilities | Hearing Impairment |
| Description | FM Receiver suitable for use with any hearing aid featuring a T-Coil that picks up the signal from transmitters. |
| Product Information Page | Link |



| | |
|---|---|
| Name of Item | Classic Inspiro FM Transmitter with iLapel mic |
| Main target users' type of disabilities | Hearing Impairment |
| Description | <p>FM Transmitter to be carried by teachers / speakers that the mic can directly pick up the signal from speaker and transmit to the FM receiver.</p> <p>(Needs to be used with a FM receiver. The RogerMyLink is available for loan and can be borrowed as one set.)</p> |
| Product Information Page | Link |



| | |
|---|--|
| Name of Item | Roger TableMic |
| Main target users' type of disabilities | Hearing Impairment |
| Description | <p>A microphone ideal for meetings or group conversations. The microphone would automatically pick up the voice of the person who are talking.</p> <p>(Needs to be used with a FM receiver. The RogerMyLink is available for loan and can be borrowed as one set.)</p> |
| Product Information Page | Link |



| | |
|---|---|
| Name of Item | Roger Pen |
| Main target users' type of disabilities | Hearing Impairment |
| Description | <p>A pen-shaped FM transmitter for portability. More versatile and can be used in various situations, for example in a social gathering whether each person passes around the pen to speak so that the student with HI can pick up each person's voice clearly through the FM Receiver.</p> <p>Supports Bluetooth connectivity and can be used for phone calls as well.</p> <p>(Note: Needs to be used with a FM receiver. The RogerMyLink is available for loan and can be borrowed as one set.)</p> |
| Product Information Page | Link |



| | |
|---|---|
| Name of Item | Zoomax Snow 8 HD Plus Handheld video magnifier (With text-to-speech function) |
| Main target users' type of disabilities | Visual Impairment |
| Description | A handheld magnifier for printed materials with image storing (up to 1000 images) and touchscreen controls. Saved images can be transferred to PCs with USB connection. |
| Product Information Page | Link |



| | |
|---|---|
| Name of Item | Recording Pens |
| Main target users' type of disabilities | Hearing Impairment / Specific Learning Difficulties |
| Description | Recording pens to record lectures |
| Product Information Page | Please click on the product names to their respective product information page. |



| | |
|---|--|
| Name of Item | C922 Pro HD Stream Webcam |
| Main target users' type of disabilities | General students |
| Description | Web camera for online classes / assessment |
| Product Information Page | Link |



| | |
|---|--|
| Name of Item | Stair Climber |
| Model Number | SUNWA TRE-52 |
| Length | 1258mm |
| Width | 757mm |
| Height | 1290mm |
| Max. capacity | 200kg |
| Main target users' type of disabilities | Physical Impairment (Wheelchair users) |
| Details | Specification |
| Notes | Please contact SEN Support in advance if you wish to use the stairclimber. |



| | |
|---|---------------------------------|
| Name of Item | Voice Amplifier Speaker |
| Model Number | CN-SBM200 |
| Main target users' type of disabilities | Students with speech impairment |
| Details | Specification |
| Notes | |